

REF: HBA -TUR -25 25 TEST

Intended Use

Hemoglobin A1c reagent is intended for Quantitative turbidimetric determination of HbA1c in human blood .

Background

The glycemic control in diabetes mellitus is mainly by the determination of glucose, but also through quantitative determination of hemoglobin A1c in human blood. HbA1c is an indication for the actual glucose levels over the preceding 3 months. It was shown that HbA1c in diabetic subjects can be elevated 2-3 fold over normal and on other hand approaches normal values when they are under metabolic control.

Assay Principle

This method utilizes the interaction of antigen and antibody to determine th HbA1c in whole EDTA blood. HbA1c in test samples is absorbed onto the surface of latex particles, whiche react with Anti-HbA1c (antigen-antibody reaction)and gives agglutination. The amount of agglutination is measured as absorbance. The HbA1c value is obtained from a calibration curve.

Reagent

Reagent1 (R1)
Latex.
Sodium azide (0.95 g/L).

Reagent2 (R2)
Anti-human hemoglobin A1c mouse monoclonal antibody.
Stabilizers.

Materials required but not provided with the kit

1- Standard set
HbA1c concentration is stated on the vials labels.

2-Controls

Reagent Preparation, Storage and Stability

HbA1c reagents are stable up to the expiry date labeled on the bottles when stored at 2 - 8°C and contaminations are prevented during their use.Once opened the reagents are stable for 1 month if stored tightly closed at 2 - 8 °C after use.

Specimen Collection and Preparation

Fresh EDTA blood.

Hemolysate procedure

To determine HbA1c, a hemolysate must be prepared for each sample as follow:

1. Dispense 2 ml hemolysis reagent into a test tube.
2. Place 20 µl of well mixed whole EDTA blood (Samples, Standards and Controls) into the test tube and mix.
3. Allow to rest 5 minutes or until complete lysis is evident.

Stability of the hemosylate: 72 hours at 2 - 8°C.

Procedure

Wavelength 650 nm
Temperature 37 °C
Cuvette 1cm light path
Zero adjustment distilled water

Solve and lyse standard/control

	Standard	Sample
Reagent (R1)	375 µl	375 µl
Standard	5 µl	-----
Sample	-----	5 µl

Mix, and incubate for 2 minutes, then add

Reagent (R2)	75 µl	75 µl
--------------	-------	-------

Mix and read absorbance (A1) immediately, then after 5 minutes read absorbance (A2).

Adaptation sheets for several automatic analyzers are available upon request.

Calculation

Generate a reference curve using HbA1c standard set. Determine Δ absorbance of the sample and each standard as following:
Δ absorbance of sample = (A2 - A1) sample
Δ absorbance of each standard = (A2 - A1) for each Standard
Plot the calibration curve and obtain the result.

Expected Values

Non-diabetics < 6 %
Therapeutic diabetics < 7 %
Each laboratory should establish its own reference range.

Linearity









Up to 15 %.
specimens showing higher concentration should be diluted 1/5 using physiological saline and repeat the assay.




Dynamic Range

0 - 15 %.

References

1. Bates, H.M., Lab. Mang., Vol 16 (Jan. 1978)
2. Gonen, B., and Rubenstein, A.H., Diabetologia 15, 1 (1978).
3. Trivelli, L.A., Ranney, H.M., and Lai, H.T., New eng. J. Med. 284, 353 (1971).

	Consult Instruction for Use
	Caution Consult Accompanying Documents
	In Vitro Diagnostic Medical Device
	Temperature Limitation
	Manufacturer
	Authorized Representative In The European Community
	Catalogue Number
	Batch Code
	Use By

 N.S BIOTEC MEDICAL EQUIPMENT 66 Port Said St., Camp Shezar Alexandria - Egypt Tele: 002 03 592 0902 Fax : 002 03 592 0908 Website : www.nsbiotec.com E- mail : info@nsbiotec.com	  CMC Medical Devices & Drugs S.L. C/ Horacio Lengo, 18. 29006. Málaga, Spain
---	--