



Ultra-Pure Water System



NS. BIOTEC
www.nsbiotec.com



Ultra-Pure Water System

- High desalinization rate: 97% - 99%
- Convenient to replace consumables.
- Safe and reliable
- On-line monitoring to ensure good quality of water output

Model	RO 100
Pretreatment cartridge	<input type="checkbox"/> sediment filter 10 inch <input type="checkbox"/> granular carbon filter 10 inch <input type="checkbox"/> bacteria filter 10 inch <input type="checkbox"/> A2 Solid carbon filter 10 inch
Membrane capacity	100 gallon/day
No. of pumps	1 pump
ultraviolet radiation unit	Less than 1 gallon/min
ion removal filter the risen tank	4*10
Reservoir tank	3 gallons
System display and indication parameters	<input type="checkbox"/> guage for inlet pressure <input type="checkbox"/> guage for outlet pressure <input type="checkbox"/> T.D.S METER
Material	the cabinet made of treated iron sheet and coated with electrostatic paints (Stainless steel optional)
Outer Dimensions (W x D x H) mm	550 x 450 x 600
Power supply	220 - 240 VAC, 50Hz



Some of the illustrations in this publication include the special accessories; we reserve the right to make changes in technical specification. Dimension without prior notice.

Headquarters : 66 Port Said St., Camp Shezar. Alexandria, Egypt
 Factory: 1 Elmadina st, Indstrial area, Elbetash, Alexandria, Egypt
 E-mail: info@nsbiotec.com export@nsbiotec.com

Tel : 002 03 5920902 Fax: 002 03 5920908

www.nsbiotec.com